CHECKPOINT®

HEAD & NECK

Model: 9394

Patented (AU, CN, DE, ES, FR, UK, HK, IT, US) -
www.checkpointsurgical.com/patents
© (AU, CA, CN, HK, IT, JP, DE, ES, FR, UK, US)

**Head & Neck Surgical Nerve Stimulator/Locator
INSTRUCTIONS FOR USE**

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CHECKPOINT®

HEAD & NECK

Model: 9394
Head & Neck Surgical Nerve Stimulator/Locator
INSTRUCTIONS FOR USE

Patented (AU, CN, DE, ES, FR, UK, HK, IT, US) -
www.checkpointsurgical.com/patents
 © (AU, CA, CN, HK, IT, JP, DE, ES, FR, UK, US)

INDICATIONS FOR USE:

The CHECKPOINT® Head & Neck Stimulator is a single-use, sterile device intended to provide electrical stimulation of exposed motor nerves or muscle tissue to locate and identify nerves and to test nerve and muscle excitability.

CONTRAINDICATION:

- ⊘ Do NOT use this Stimulator when paralyzing anesthetic agents are in effect, as an absent or inconsistent response to stimulation may result in inaccurate assessment of nerve and muscle function.

WARNINGS:

- ⊘ Do NOT apply stimulation across the patient's chest because the introduction of electrical current into the chest may cause rhythm disturbances to the patient's heart.
- ⊘ Do NOT leave Stimulator unattended in the surgical field.
- ⊘ Do NOT use Stimulator in the presence of flammable anesthetics.
- ⊘ Do NOT use this Stimulator in patients with implanted electronic devices. Stimulus current near such devices may cause them to malfunction.
- ⊘ Do NOT use in cerebrospinal fluid, brain, meninges, or spinal cord.
- ⊘ Do NOT use the Stimulator while delivering electrocautery. Electrocautery may interfere with or damage the Checkpoint Head & Neck Stimulator or cause injury.
- ⚠ Use of a tourniquet may reduce nerve and muscle excitability distal to the tourniquet. If a decrease in nerve or muscle response to stimulation is observed, it may be necessary to take down the tourniquet and allow time for re-perfusion before testing with Checkpoint.
- ⊘ Do NOT attempt to reuse with another patient. Serious patient cross-infection or pyrogenic reaction could occur.

**DO NOT APPLY
ACROSS CHEST**

PRECAUTIONS:

- Single use only. Once turned ON, the Stimulator will not turn OFF and will remain operational for approximately four (4) hours.
- Use the Stimulator at 0.5mA initially and increase the amplitude as needed.
- ⊘ Do NOT bend the probe or apply excessive pressure where the probe enters the nose cone. Do NOT use the Checkpoint Head & Neck Stimulator as a dissecting tool or for any other purpose not specified in these Instructions for Use.
- ⚠ The needle return electrode may cause stimulation. Place the needle return electrode in an area away from excitable tissue and away from the stimulation area to avoid confusion.
- ⊘ Do NOT submerge or saturate the Checkpoint Head & Neck Stimulator with fluid.
- ⚠ The Stimulator may interfere with patient monitoring equipment in the operating room while stimulating.
- ⊘ Portable and mobile Radio Frequency (RF) communications equipment can interfere with the Stimulator. Do not use such equipment while delivering stimulation.
- The Stimulator is not likely to cause any interference in nearby electronic equipment through Electromagnetic Interference (EMI) (complies with CISPR 11, Type A). For additional information on Electromagnetic Interference:
<http://www.checkpointsurgical.com/510k/>
- .
- ⚠ Be aware that using the device in the vicinity of certain nerves can have systemic effects (e.g. bradycardia, hypotension, laryngospasm).

The Checkpoint Head & Neck Stimulator is an adjunctive tool for locating nerves and testing nerve and muscle excitability. It is not intended to replace good surgical techniques in locating and dissecting tissue. The surgeon is ultimately responsible for identifying exposed motor nerves and tissue.


INSTRUCTIONS FOR USE:

1. Prepare the Checkpoint Head & Neck for use

- Remove the sterile pouch from the package.
- Inspect sterile pouch for any damage.
- Do NOT use if sterile pouch or product appears damaged.
- Remove product from sterile pouch.
- Remove protective sleeve from probe.

- Remove tape from return electrode cable.
- Remove cover from the needle return electrode.

Do NOT turn on until ready to use. Once turned ON, the device will not turn OFF and will remain operational for approximately four (4) hours.

2. Place the needle return electrode in subcutaneous tissue:

- In the sterile field but away from the area to be stimulated. If stimulation occurs at the needle site, move the needle.
- Do NOT cross midline of the body.
- Do NOT place in muscle or near other excitable tissue.

The needle return electrode may cause stimulation. Place the needle return electrode in an area away from excitable tissue and away from the stimulation area to avoid confusion.

3. Prepare to deliver stimulation

Probe	Nose Cone	Stimulus Amplitude Switch
↓	↓	↓

- **NOTE: The Stimulator must be turned ON BEFORE the probe makes contact with tissue.**
- Turn ON the Checkpoint by moving the stimulus amplitude switch from 0 to 0.5mA.
- The visual indicator will turn Flashing RED, indicating that the Stimulator is ON but the circuit is not complete.
- Increasing amplitude increases the strength of stimulation.
- Amplitude switch provides three (3) settings 0.5mA @ 50 μsec, 1.0mA @ 100 μsec, 2.0mA @ 150 μsec.
- If a nerve is thought to be in close proximity and no response is obtained at the lowest amplitude settings, the surgeon may increase the amplitude to 1.0mA then to 2.0mA, if needed.
- Pulse duration is preset at each of the amplitude settings.
- The stimulator delivers a continuous train of stimulation, which will result in a sustained response as long as the stimulus is applied.

4. Deliver stimulation

- Place probe in contact with tissue.
- The visual indicator will FLASH YELLOW indicating that stimulation is being delivered.

NOTE: If a questionable response to stimulation is obtained from exposed motor nerves or muscle tissue under test, move the probe to a known nerve or muscle and observe the results to confirm that the Stimulator is operating correctly. Check connections, check controls, and move around the area. Verify that any paralytic agents used are no longer in effect.

Additional notes:

- A FLASHING RED indicator means that stimulation has been requested, but adequate stimulus current is NOT being delivered because of poor connection of the probe or the needle return electrode to the patient tissue.
 - Check needle return electrode contact and position.
 - Check probe contact and position.
 - Move probe to a position with a known response.
- Moving the stimulus amplitude switch to the OFF/Standby position will stop stimulation and the visual indicator will be illuminated YELLOW continuously.

End of use / disposal:

- A low (depleted) battery will cause the Stimulator to automatically turn off and the visual indicator will not be illuminated. No further use of the Stimulator will be possible.
- At end of use, move the stimulus amplitude switch to the OFF/Standby position.
- Cut off and dispose of needle electrode in appropriate sharps / biohazard container.
- Product is powered by two size N alkaline cells in series permanently sealed within the device. Use appropriate disposal methods per hospital guidelines.


INDICATOR LIGHT STATUS:

LIGHT	STIMULUS STATUS
Solid yellow	The Stimulator has been turned on but the amplitude switch is now in the Off/standby position. Stimulation is NOT being delivered.
Flashing yellow	Stimulation is being delivered (NOTE: the flashing rate does NOT correspond to the stimulus frequency or intensity).
Flashing red	Stimulation has been requested, but adequate stimulus current is NOT being delivered because of poor connection of the probe or the needle return electrode to the patient tissue.

WARRANTY and LIMITATIONS:

Checkpoint Surgical, INC. (CPS) warrants that this product has been manufactured, packaged, and tested with reasonable care and will be free from defects in workmanship and materials. CPS further warrants that the product will remain sterile for a period described on the product's label, provided the original packaging remains intact. This product is for single-use only and is not intended or designed for reuse. This warranty shall not apply to product that has been re-sterilized, repaired, altered, or modified in any way, or to products that have been improperly stored or operated. CPS will not be liable for any incidental, special or consequential loss, damage, or expense resulting, directly or indirectly, from the use of the product. The sole obligation of CPS shall be to refund or replace, at its option, any device that CPS determines was defective at time of shipment if notice thereof is received before expiration date described on such product label. Buyer assumes all liability, whether based upon warranty, contract, negligence, or otherwise, for damage resulting from the handling, possession, use or misuse (including reuse) of this product. Because CPS has no control over the operation, inspection, maintenance, or use of its products after sale and has no control over selection of patients, THIS WARRANTY IS EXPRESSLY IN LIEU OF ANY OTHER EXPRESS OR IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE, AND OF ANY OTHER OBLIGATION ON THE PART OF THE SELLER. The remedies set forth in the Warranty and Limitations shall be the exclusive remedy available to any person. No agent, employee, or representative of CPS has any authority to change any of the foregoing or assume or bind CPS to any additional liability or responsibility in connection with this warranty. Buyer's use of this product shall be deemed acceptance of the terms and conditions of this Warranty and Limitations.

CHECKPOINT (P/N 9394) SPECIFICATIONS:

Stimulus Frequency	16Hz (+/- 1Hz)
Stimulus Amplitude	User Selectable: 0.5mA (+/- 0.05mA)@ 50±4 μsec 1mA (+/- 0.1mA)@ 100±4 μsec 2mA (+/- 0.2mA)@ 150±4 μsec
Stimulus Waveform	Biphasic stimulus waveform with controlled current during the cathodic (leading) phase and no net DC current. Specified operation is ensured into any patient circuit load resistance up to 1,500Ω. 
Operating Life	Approximately four (4) hours of operation
Power Source	Two size N alkaline cells (permanently sealed within device)
Operating Temperature	+16°C to +26°C (61°F to 79°F)
Relative Humidity	30% to 60%RH
Atmospheric Pressure	88kPa to 106kPa
Storage Temperature	-20°C to +55°C

- Product has no user serviceable or repairable parts.
 - Sterilized with Ethylene Oxide.
 - Protection against water ingress: Ordinary Equipment (IPX0).
-  Do NOT resterilize.

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Rx only

EN This system does not contain natural rubber (including natural rubber latex, dry natural rubber, synthetic latex and synthetic rubber that contains natural rubber in its formulation).

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EN Caution: U.S. Federal Law restricts this device to sale by or on the order of a physician or properly licensed practitioner.

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*+B1309394u/
\$\$3YYMMDDL*LL*.

EN UDI barcode with Checkpoint LIC, part number, unit level (1: sterilized pouch, 2: single unit sales carton), expiration date and lot number.

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
CE
2797

Manufactured for:
Checkpoint Surgical, Inc.
6050 Oak Tree Blvd
Independence, OH 44131 USA

www.checkpointsurgical.com

Winckels Medical Devices Expertise
WMDE B.V.
Bergerweg 18
6085 AT Horn
The Netherlands
office@wmdegroup.com

 **WARNING: Cancer
and Reproductive Harm -
www.P65Warnings.ca.gov**

© Checkpoint Surgical, Inc. Jan 2023